




Wall Mount Cup Sinks



SYMBOL	OUTLET	A	B	C	D	E	F	G	H	J	THREAD	K
W39026-150	1.5 NPS	5.88"	2.63"	4.63"	7"	3.88"	5.13"	6.88"	4.63"	2"	1.5"	.25"


## Wall Mount Cup Sink With Backnut (Included)



SYMBOL	OUTLET	A	B	C	D	E	F	G		K
W39027-150	1.5 NPS	5.88"	2.63"	4.63"	7"	3.88"	5.125"	6.88"	1.5"	.25"

BACKNUT = W38953-166

## Wall Mount Cup Sink



SYMBOL	OUTLET	A	B	C	D	E	F	G	H	J	THREAD	LIP
W39034-150	1.5 NPS	3.18"	3.18"	3"	4"	5.18"	4.5"	4.75"	4.75"	2"	1.5"	.25"
W39034-151	1.5 NPS	3.18"	3.18"	3"	4"	5.18"	4.5"	4.75"	4.75"	2"	1.5"	.25"
W39035-150	1.5 NPS	5.50"	6.56"	4"	6.5"	19"	4.5"	6"	4.75"		1.5"	.25"
W39035-151	1.5 NPS	5.50"	6.56"	4"	6.5"	19"	4.5"	6"	4.75"	2"	1.5"	.25"