


## Center Drain Sinks



This series of polyolefin sinks allows uniform drainage from all locations on the sink. Steel rims available for undercounter mount. Stainless steel rims available for flush mount. Standard color is non-glare black. Sinks available in other colors on special order. Contact your sales professional for details.

SYMBOL	A	B	C	D	E	F	G	H	OUTLET
W31415-000	11"	11"	7"	13.5"	13.5"	5.5"	9.75"	9.75"	1.5"
W31418-000	12"	8"	8"	14.5"	10.5"	4"	3.5"	11"	1.5"
W31424-000	13.75"	4.5"	5.5"	15.375"	6"	2.25"	6.75"	7.1875"	1.5"
W31440-000	16"	6"	6"	18"	8"	2.75"	7.75"	8.5"	1.5"
W31449-000	16"	6"	6"	18.5"	8.5"	2.75"	7.75"	8.5"	1.5"
W31466-000	18"	15"	12"	20.5"	17.5"	7.5"	9"	15"	1.5"
W31473-000	21"	18"	12"	23.5"	20.5"	9"	10.5"	15"	1.5"
W31478-000	23"	17"	11.5"	25.5"	19.5"	8.75"	11.5"	14"	1.5"
W31481-000	24"	20"	10"	26.5"	22.5"	10"	12"	13"	1.5"
W31483-000	24"	20"	12"	26.5"	22.5"	10"	12"	15"	1.5"